Recrystallization of Salicylic acid

Recrystallization

Solid organic cpd.s when isolated from organic reaction are impure; they are contaminated with small amounts of other cpd.s produced along with the desired product.

The purification of impure crystalline compound is usually done by Recrystallization from a suitable solvent or a mixture of solvents.


Purification of solids by recrystallization is based upon differences in their solubility in a given solvent or a mix. of solvents.

Solid organic cpd.s when isolated from organic reaction are impure; they are contaminated with small amounts of other cpd.s produced along with the desired product.

Choosing a solvent for Recrystallization:

The ideal solvent should:

- 1- Chemically inert toward the solute.
- 2- It should dissolve the solute to be purified readily at or near it's boiling point, but sparingly at the lab. temp. or below (0 25 °C).
- 3- It should dissolve the impurities readily or not at all.
- 4- not be flammable, low cost, and of low toxicity.




Simple Recrystallization process consist of:

- 1- Dissolving the impure substance in some suitable solvent at or near the boiling point.
- 2- Filtering the hot solution from the particles of insoluble material & dust, (Hot filtration).
- 3- Allowing the hot solution to cool thus causing the dissolved substance to crystallize out.
- 4- Separating the crystals from the supernatant soln.

How could we choose a good solvent:

Practically,

Take 0.1g of a pure sample of cpd. to be purified & try to dissolve it in 1ml of solvent,





Results of solubility tests for cpd. (A) are shown in table (g/ml).

Solvent	Water	Ethanol	Diethyl ether
Cold	20	3	5
Hot	30	25	5

Which solvent will you choose to recrystallize cpd. A?

Using Charcoal:

Samples to be purified may contain soluble colored impurities that may cause the soln. & the crystals to be colored.

Up on Recrystallization these impurities dissolves in the boiling solvent & adsorbs on the crystals produced up on cooling yielding a colored product.

Activated charcoal composed of fine carbon particles with a large active surface area on which the colored impurities will be adsorbed.

Charcoal is added to the hot soln. before boiling & the soln. is kept hot at or near the b.p. for about 3 - 5 min. with shaking to wet the charcoal, the solution is then filtered through a fluted filter paper.

V

Q / Charcoal is Not used for recrystallization of phenolic cpd.s
Ans: Because, They contain ferric ions (Fe) that upon heating
the solution for some times it can react with the phenolic -OH
group forming red - violet colored complexes thus impairing the
purification process.


Recrystallization using mixed solvents:

It is applied when our cpd. is readily soluble in a solvent at room temp. & insoluble in other solvent, The 2 solvents must be miscible with each other as Alcohol & water, ether & pentane glacial acetic acid & water.

Name of Experiment: Purification of the synthesized Aim of experiment: Purification of the synthesized Salicylic Acid.

Procedure:

- 1- Put 1 g of impure Salicylic acid sample in a beaker.
- 2- Try to dissolve it in a minimum amount of hot Continuously water with heating.
- 3- Filter the solution while it is hot.
- 4- Cool the filtrate, then Salicylic acid will crystallize.
- 5- Filter again, (Cold filtration).
- 6- Collect the crystals of S.A. on the filter paper and dry them on oven.



Calculations:

```
Percent pure substance % = 

Sample wt. x 100
```

Percent impure substance % = 100 - Percent pure substance